

COSMOLUX[®] VHR[®]


COSMOLUX[®] VHR[®]



| | Product Number | Description | Wattage (W) | Base | UVA Flux (W) | UVB Flux (W) | Total UV (W) | Rated Hours | Typical Exposure Time |
|------------------------|----------------|-------------|----------------|---------------|--------------|--------------|--------------|-------------|-----------------------|
| CosmoLux VHR XT | 16394 | FR71 | 160-200 | Bi-Pin | 41.8 | 1.3 | 43.1 | 800+ | 7-12 Mins. |
| CosmoLux VHR TT | 16239 | FR74 | 160-200 | RDC | 38.0 | 1.2 | 39.2 | 800+ | 7-12 Mins. |
| CosmoLux VHR TT | 16107 | FR73 | 160 | RDC | 38.0 | 1.2 | 39.2 | 800+ | 7-12 Mins. |
| CosmoLux VHR TT | 16169 | FR71 | 160 | Bi-Pin | 38.0 | 1.2 | 39.2 | 800+ | 7-12 Mins. |
| CosmoLux VHR TT | 16112 | FR59 | 140 | Bi-Pin | 31.0 | 1.0 | 32.0 | 800+ | 7-12 Mins. |
| CosmoLux VHR | 16220 | FR2M | 180-200 | Bi-Pin | 41.5 | 1.2 | 42.7 | 800+ | 9-15 Mins. |
| CosmoLux VHR | 16180 | FR71 | 160 | Bi-Pin | 38.0 | 1.1 | 39.1 | 800+ | 9-15 Mins. |

VHR XT.
The most robust standard spectrum lamp in the industry.

**FR71 Cosmolux VHR-XT 160-200W
 Spectral Power Distribution**



Contraindication: This product is contraindicated for use on persons under the age of 18 years.
 Contraindication: This product must not be used if skin lesions or open wounds are present.
 Warning: This product should not be used on individuals who have had skin cancer or have a family history of skin cancer.
 Warning: Persons repeatedly exposed to UV radiation should be regularly evaluated for skin cancer.