

COSMOLUX[®] LYNX SERIES

COSMOLUX[®] LYNX[®]

	Product Number	Description	Wattage (W)	Base	UVA Flux (W)	UVB Flux (W)	Total UV (W)	Rated Hours	Typical Exposure Time
CosmoLux Lynx Series S1	26500	FR71	120	Bi-Pin	28.0	1.0	29.0	800+	10-12 Mins.
CosmoLux Lynx Series S2	26501	FR71	120	Bi-Pin	28.0	1.3	29.3	800+	10-12 Mins.
CosmoLux Lynx Series S3	26502	FR2M	200	Bi-Pin	41.5	2.2	43.7	800+	10-12 Mins.
CosmoLux Lynx Series K1	26503	FR75	180	Bi-Pin	39.8	1.4	41.2	800+	10-12 Mins.
CosmoLux Lynx Series K2	26504	FR71	200	Bi-Pin	45.5	2.0	47.5	800+	10-12 Mins.
CosmoLux Lynx Series S4	26506	FR71	160	Bi-Pin	42.6	1.0	43.7	800+	10-12 Mins.
CosmoLux Lynx Series U1	26507	FR79	180	Bi-Pin	41.7	1.1	42.7	800+	10-12 Mins.

COSMOLUX[®]
LYNX

COSMOLUX[®]
LYNX

COSMOLUX[®]
LYNX

COSMOLUX[®]
LYNX

Contraindication: This product is contraindicated for use on persons under the age of 18 years.
 Contraindication: This product must not be used if skin lesions or open wounds are present.
 Warning: This product should not be used on individuals who have had skin cancer or have a family history of skin cancer.
 Warning: Persons repeatedly exposed to UV radiation should be regularly evaluated for skin cancer.